



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

## Wood Vinyl 240 Capacity Punch Card Medication Cart, Standard Package #WV240PC-OAK

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.



### Specifications

#### Construction:

- Vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate

#### Features:

- Light weight vinyl coated aluminum
- Single column medication cart with 240 punch card capacity
- Key locking cabinet with locking narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- 5" ball-bearing premium rubber casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture

#### Skirt Option:

- #WV-2WSKIRT-OAK – Skirt option to cover wheels

#### Drawer Configuratin:

- 1 ea – 3" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers

#### Assembled Cart Dimensions:

- 44.5" H x 19" W x 25" D
- Writing Surface: 44"
- Cart weight: 102 lbs
- Approx shipping weight: 122 lbs



Skirt Option



Cherry Mahogany (CM)

WV240PC-CM



Buffalo Oak (OAK)

WV240PC-OAK

#### Color Options:

The WV240PC-OAK is shown above in buffalo oak finish with a grey granite finish on the top. Also available in cherry mahogany #WV240PC-CM  
 Additional vinyl color options available with minimum cabinet purchase and extended lead time.