

ADA305AF



32" x 19.75" x 22.5" (H x W x D)

Built-in undercounter -25°C ADA compliant commercially-approved all-freezer in white with lock, digital controls, interior baskets, hospital cord with 'green dot' plug, factory installed access port, and manual defrost operation

Highlights:

Medical grade design to function at lower temperatures (-25°C) than standard residential units

Just 32" high for an easy fit under lower ADA compliant counters

Slide-out storage baskets with plastic fronts help to protect stored contents from warmer ambient temperature

Product Features:

ADA compliant design	32" height allows easy installation under lower ADA compliant counters
-25°C capable	Low temperature operation ideal for vaccines, pharmaceuticals, and other sensitive materials
Commercially approved	Meets ETL-S to NSF-7 standards for use in commercial settings
Built-in capable	Front-breathing design lets you make the best use of space by installing your appliance under the counter
Factory-installed lock	Located towards the top of the door for easy access and to provide added security of unit (2 keys included)
Access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Adjustable digital thermostat	Electronic controls located in the front kickplate for easy & institutive temperature management, with a readout in Fahrenheit or Celsius
Basket shelves	Three slide-out baskets offer convenient storage, with plastic fronts to help protect contents from warmer air when the door is opened

3" thick door	Improved insulation for an efficient performance; NOTE: due to the oversized design, the door does open in its own radius. If installing this freezer against a wall, please see our drawings to ensure adequate space
---------------	--

ADA305AF Specifications:

Overview	
Height of Cabinet	32.0" (81 cm)
Height to Hinge Cap	32.0" (81 cm)
Width	19.75" (50 cm)
Width with Door Open	20.25" (51 cm)
Depth	22.5" (57 cm)
Depth with Handle	24.38" (62 cm)
Depth with door at 90°	40.25" (102 cm)
Capacity	2.47 cu.ft. (70 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	65.0 lbs. (29 kg)
Shipping Weight	70.0 lbs. (32 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
Freezer	
Door Swing	RHD
Reversible	Yes
Interior Height	23.75" (60 cm)
Interior Width	14.5" (37 cm)
Interior Depth	17.0" (43 cm)
Shelf Type	Drawers
Shelf Qty	3
Thermostat Type	Digital
Fan Type	Exterior
Refrigerant Type	R600a
Refrigerant Amount	1.34 oz.
High Side PSI	180.0
Low Side PSI	80.0
Level Legs	2
Interior Drawers	Yes
Drawers Qty	3
Compressor Step Height	5.75" (15 cm)

Compressor Step Width	14.5" (37 cm)
Compressor Step Depth	6.0" (15 cm)
Temperature Range	-25° to -15°C

