

## Open Base Power Table with Adjustable Backrest, Footrest and Stirrups

### MODEL 81360

- Pneumatic backrest is adjustable from either side and raises 80 degrees from flat
- Footrest is adjustable from either side and raises 80 degrees from flat
- Retractable, 2-position, length adjustable stirrups
- Easy-clean, ABS plastic enclosed lift
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 2" ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 500 lbs. (227 kg) load capacity under normal use



### 81360 SPECIFICATIONS

Overall	76" L x 27" W x 18–35" H
Backrest	38" L x 27" W
Footrest	13" L x 19" W

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



### Upholstery Colors



## OPTIONS

### 080 Safety Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 safety straps per set (92" x 2")

080



### 082 Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



### 035 Armboard

- Padded, ratchet style, rotating armboard
- Rotates 180°
- Interchangeable on right or left side
- 7" W x 17.5" L

035



### 012 Headrest

- Padded and length adjustable
- Extends the table's length
- 11" tapers to 9" x 10" H

### 084 Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

012



084



087



### 087 Locking Casters

- Heavy duty, rubber wheel, easy-roll casters


### 030 Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030



040



### 040 Paper Cutter

- Clear plastic with elastic straps and buckle

### 034 Paper Dispenser/Cutter Combo

- Paper dispenser and paper cutter ordered as one option