



M-series

MDS3027B05

M-SERIES EMERGENCY CART with Five Drawers and a Breakaway Lock

Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel break-away lock bar
- Pack of 100 break-away seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- 27" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

Engineered &
Manufactured
in the U.S.



Color shown is Red. More than
20 additional colors are available.



Medium Height, Standard Width, Five Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

Drawer Configuration:

- Two 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.9 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

Assembled Cart Dimensions:

- 40.5" H x 36.75" W x 22" D (102.9 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 145 lbs (65.8 kg)
- Approximate shipping weight: 184 lbs (83.5 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty:

12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/acute-care-carts/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300