



# M-series MDS3030B06+MD30-EMG3

## M-SERIES TALL EMERGENCY CART with MD30-EMG3 Emergency Package

### MD30-EMG3 Accessories Include:



Oxygen tank holder (O2HLDR)



Large Utility Basket (BASKETDM)



(2) Utility hooks (UHOOKS2)



Engineered & Manufactured in the U.S.



Toprail, matching cart color (MD30-TOPRAIL)



Cardiac board and brackets (MD-CARDBRD)



Chart Holder (CHARTHLDRDM)



(2 Sets) Standard 6" Drawer Dividers (MD30-DIV6-S)



(1 Set) Standard 3" Drawer Dividers (MD30-DIV3-S)



(1 Set) Standard 3" Drawer Tray Dividers (MD30-TRAYDIV3-S)




IV pole (MD-IV2)



Defibrillator Shelf (MD-DEFIB)



Locking Sharps Container (5qt) & Glove Box Holder (LKSHRPHOLD5Q-GLV)



Suction shelf (SUCTIONSHLF)

**Color shown is Red. More than 20 additional colors are available.**



#### Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

#### Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D) cart only,
- 71.81" H x 48.06" W x 22" D (182.4 cm H x 122.1 cm W x 55.9 cm D) w/accessories
- Cart weight: 213 lbs (96.6 kg)
- Approximate shipping weight: 252 lbs (114.3 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

**Toll Free: (800) 433-4064 • Ph: (719) 637-0300**  
**medicalsales@harloff.com • www.harloff.com**

**HARLOFF**  
INNOVATING STORAGE FOR CARE