



M3DS1824K04

M-SERIES SHORT MINI CART

Four Drawers with a Standard Key Lock

Features:

- · Sturdy, uni-body steel construction
- · Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- · Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- · Antimicrobial plastic top, removable
- · 24" of vertical drawer space
- · Cart is delivered assembled
- · Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- · Latex free

Colors:

Shown in Yellow, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 33.375" (84.8 cm)
- Cart dims: 34.5" H x 18" W x 18" D (87.6 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 83 lbs (37.6 kg)
- Approx shipping weight: 107 lbs (48.5 kg)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm H x 61 cm W x 61 cm D)

Drawer Configuration:

One 3" drawer (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Short Height, Mini Width, Four Drawer Infection Control Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.





View similar products online at

https://www.harloff.com/product-category/m-series-medical-carts/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300