



MPA1830B06



A-SERIES MINI WIDTH CART

with Six Drawers and a Breakaway Lock

Features:

- Extremely light weight cart with critical grade aluminum construction, ideal for carts that are moved consistently throughout the day
- · Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- · Antimicrobial plastic top, easy to remove/replace
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling Tente® casters one steering caster and two braking casters
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at https://www.harloff.com/locking-options/
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- · Cart is delivered assembled
- · 30" of vertical drawer space

Engineered & Manufactured in the U.S.

Tall Height, Mini Width, Six Drawer Cart with a Breakaway Lock

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.9 cm)
- 12.5" W x 17" D (31.8 cm x 43.2 cm)

Assembled Cart Dimensions:

- 43.75" H x 24.75" W x 22" D (111.1 cm H x 62.9 cm W x 55.9 cm D)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm x 61 cm x 61 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



https://www.harloff.com/product-category/emergency-crash-carts/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300



