

aluminum





A-SERIES STANDARD WIDTH CART

with Five Drawers and a Breakaway Lock

HARLOFF

Features:

- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest guality in the industry
- Heat-sealed durable scratch resistant
 powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Antimicrobial plastic top, removable and easily replaceable
- Ergonomic side push handles can be mounted on the left or right
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling Tente® casters -1 steeringcaster and 2 braking casters
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at <u>https://www.harloff.com/locking-options/</u>
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- · Cart is delivered assembled
- 24" of vertical drawer space

Engineered & Manufactured in the U.S.

Short, Standard Width, Five Drawer Crash Cart with a Breakaway Lock

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices. Color shown with a Red body and drawers. More than 20 additional colors are available.

Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Three 3" drawers (7.6 cm), one 6 " drawer (15.2 cm) and one 9" drawer (22.8 cm)

Assembled Cart Dimensions:

- 37.25" H x 36.75" W x 22" D (94.6 cm H x 93.3 cm W x 55.9 cm D)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts.



View similar products online at

https://www.harloff.com/product-category/emergency-crash-carts/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300