

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Universal Line Wide Medication Cart with Key Lock, Standard Package (MEDICINE ON TIME® Compatible) #AL8W10K4MOT

Harloff's Universal Line of medication carts were designed to accommodate all types of medication packaging. Its modern appearance incorporates industry-leading standard features and offers multiple locking options.



Seventeen available color choices:

(may be viewed at www.harloff.com)

Beige, Burgundy, Light Blue, Navy, Forest Green, Light Gray, Mauve, Red, Teal, Yellow, White, Hammer Black, Hammer Red, Hammertone Blue, Hammertone Gray, Hammertone Green, and Sand

Specifications

Construction:

All Aluminum construction

Features:

- Extended width single column medication cart with up to 250 MOT® card capacity
- · Locking narcotics storage included
- Internal pullout shelf on right side
- Key locking cart with 4 keys per lock
- Painted steel top rail / push handle
- · Removable plastic top
- Protective pontoon bumper
- Extended caster footprint for stability when drawers are extended
- 5" ball-bearing premium casters, all full swivel, 2 with toe locking brake & 1 directional
- Ball –bearing full extension drawer slides
- · Cart is delivered fully assembled and ready to use
- Pre-drilled hole locations for standard accessories
- Latex free
- Durable, hard-baked power coat finish available in 17 colors
- Twelve year manufacturer's warranty

Drawer Configuration:

- 1 ea 3.25" multipurpose storage drawer
- 1 ea 10" drawer with one aluminum riveted divider and locking narcotics box.
- 2 ea 10 drawer with two aluminum riveted dividers.

Internal Drawer Dimensions:

• 30" W x 17" D

Assembled Cart Dimensions:

- 48.5" H x 38.25" W x 25" D
- Writing Surface: 44"
- Cart weight: 125 pounds
- · Approximate shipping weight: 160 pounds