



MPA3027E06



Color shown with

a Beige body and

Navy drawers.

More than 20

additional colors

are available.

A-SERIES STANDARD WIDTH CART


with Six Drawers & Electronic Pushbutton Lock


Features:

- · Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- · Heat-sealed durable scratch resistant powder coat finish, inside and out
- · Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- · Antimicrobial plastic top, removable and easily replaceable
- · Ergonomic side push handles can be mounted on the left or right
- · Fully customizable drawer configurations at no extra charge
- · Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Unmatched selection of locking options. types on the same cart, alarm access locks at https://www.harloff.com/locking-options/
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- · Cart is delivered assembled



Engineered & Manufactured in the U.S.







Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Four 3" drawers (7.6 cm), one 6"

drawer (15.2 cm) and one 9" drawer (22.7 cm)

Assembled Cart Dimensions:

- 40.5" H x 36.75" W x 22" D (102.9 cm H x 93.3 cm W x 55.9 cm D)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Medium Height, Standard Width, Six Drawer Cart with Keyless **Electronic Pushbutton Lock**

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



View similar products online at

https://www.harloff.com/product-category/anesthesia-carts/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300