

ARS8PV-FS30LSTACKMED2



76.5" x 23.38" x 24.38" (H x W x D)

Stacked combination of ARS8PV all-refrigerator with antimicrobial silver-ion handle and hospital grade cord with 'green dot' plug and FS30LMED2 compact manual defrost all-freezer for vaccine storage

Highlights:

Independent refrigerator and freezer storage in one slim-fitting footprint for space-challenged facilities

Intelligent microprocessor digital temperature controller for optimum temperature management

ARS8PV refrigerator has a powder-coated handle with naturally occurring silver ions to help reduce germ spread

Product Features:

| | |
|-------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Space-saving design | Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint |
| ARS8PV all-refrigerator | 8 cu.ft. automatic defrost all-refrigerator with digital controls, wire shelving, and keyed lock for added security |
| FS30LMED2 all-freezer | Compact manual defrost all-freezer with digital thermostat and wire shelving |
| Antimicrobial silver ion-coated handle (refrigerator) | Powder-coated handle with naturally occurring silver ions to help reduce germ spread |
| Easy assembly | Unit ships in three boxes for customer to install the StackRack on the bottom refrigerator/freezer and place and secure the other unit on top |
| Purpose-built design | Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC guidelines |
| +/-1°C Variation (refrigerator) | The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period |
| Adjustable temperature | Intuitive controls let you adjust the interior temperature between +2 to +8°C (refrigerator) and -20° to -15°C (freezer) |

Temperature alarms, Visual and audible alarms alert if the interior temperature rises or falls out of set range

ARS8PV-FS30LSTACKMED2 Specifications:

| Overview | |
|-----------------------------------|---------------------|
| Height of Cabinet | 76.5" (194 cm) |
| Width | 23.38" (59 cm) |
| Width with Door Open | 23.88" (61 cm) |
| Depth | 24.38" (62 cm) |
| Depth with Handle | 26.13" (66 cm) |
| Depth with door at 90° | 46.5" (118 cm) |
| Capacity | 9.8 cu.ft. (278 L) |
| Defrost Type | Automatic |
| Door | White |
| Cabinet | White |
| US Electrical Safety | ETL |
| Canadian Electrical Safety | ETL-C |
| Amps | 2.2 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Weight | 204.0 lbs. (93 kg) |
| Shipping Weight | 240.0 lbs. (109 kg) |
| Parts & Labor Warranty | 2 Years |
| Compressor Warranty | 5 Years |
| UPC | 761101077659 |
| Refrigerator-Freezer | |
| Door Swing | RHD |
| Reversible | Factory Reversible |
| Adjustable Shelves | Yes |
| Interior Light | Yes |
| Refrigerator - Interior Height | 41.25" (105 cm) |
| Refrigerator - Interior Width | 20.0" (51 cm) |
| Refrigerator - Interior Depth | 17.0" (43 cm) |
| Refrigerator - Capacity | 8.0 cu.ft |
| Refrigerator - Shelf Type | Wire |
| Refrigerator - Shelf Qty | 7 |
| Refrigerator - Temperature Range | 2 to 8°C |
| Refrigerator - Refrigerant Type | R600a |
| Refrigerator - Refrigerant Amount | 1.73 oz. |
| Freezer - Interior Height | 18.5" (47 cm) |

| | |
|------------------------------|---------------|
| Freezer - Interior Width | 13.5" (34 cm) |
| Freezer - Interior Depth | 13.0" (33 cm) |
| Freezer - Capacity | 1.8 cu.ft |
| Freezer - Shelf Type | Wire |
| Freezer - Shelf Qty | 2 |
| Thermostat Type | Digital |
| Fan Type | Interior |
| Freezer - Temperature Range | -20° to -15°C |
| Freezer - Refrigerant Type | R600a |
| Freezer - Refrigerant Amount | 0.8 oz. |

