

**FOR ENT SCOPES**

*and other small scopes*



**SureDry**  
SCOPE CABINETS

**6 ENT SCOPE CAPACITY WALL MOUNT  
DRYING CABINET, KEY LOCK**

**SCW2430DRDP**

*This new wall mount ENT scope cabinet allows for ENT scopes to be dried in the same safe, efficient manner of scopes in other parts of the facility*

- HEPA filtered Positive Pressure with two High Volume intake fans and one exhaust fan
- Customizable scope holders to fit most common scope types
- Holds up to 6 ENT Scopes
- Glass paneled doors for easy viewing of scope inventory
- Adjustable timer (10 to 180 mins) for fan operation.
- Fans activated automatically upon closing of doors
- Industry best 12 year limited manufacturer warranty
- Ideal for non-connected or removable video/light processor cord
- Cabinet dimensions: 30"H x 24"W x 12"D

*Comes in Stainless Steel or heat sealed powder coat low carbon steel, available in over 20 color options*

**Engineered & Manufactured in the U.S.**



**This optional LED-UV disinfection strip helps to keep scopes sterilized!**



**UV**

**HARLOFF**  
INNOVATING STORAGE FOR CARE

UV-C disinfection range is from 185 to 280nm. Our LED runs at 270nm, close to the optimal disinfection wavelength.

View similar products online at

<https://www.harloff.com/scope-drying-and-storage-cabinets/suredry-scope-cabinets/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**



Scopes not included

# 10 Reasons to Choose **HARLOFF** for Your Scope Cabinets

*A Smart Choice to Meet & Exceed Current SGNA, AORN and AAMI Guidelines*



1

**PureDry SCOPE CABINETS**  
*Industry Best Positive Pressure HEPA Filtered Drying System.* Designed to add Dri-Scope Aid® now or later with easy “Plug & Play” installation.

## Proudly Made in the USA

6



Harloff products are proudly **Designed, Engineered and Manufactured in the United States** at our Colorado Springs, CO facility.

## Optional LED UV-Disinfecting Lights


2

*Help to keep your scopes sterilized* with an optional UV radiation light strip as a way to disinfect your surfaces. Lights trigger in tandem with the fans after the cabinet doors are closed for scope drying.

## Easy Access to HEPA Filters

7

*Change HEPA filters easily with just thumb screws.* Save time and minimize cycle interruptions.



## Optional Wall or Cart Mounting

3

*Our cabinets fit your floor plan!* Use our wall-mounted cabinet for permanent installations or choose a mobile cart-mounted cabinet to move scopes from your reprocessing area directly to your treatment room.

## Keyed or Keyless Locks

8

*Options include* latch locks, key locks or standard electronic push button locks. View your locking options online at [www.harloff.com/locking-options/](http://www.harloff.com/locking-options/)



## Versatile Scope Holders

4

Holds major brands including **Optim, Fuji, Olympus, Pentax, Karl Storz and others.**



## Scope Cabinet Finish Options

9

Choose Stainless Steel, Heat Sealed Powder Coat in High Gloss or Textured Low Gloss, or upgrade to an **antimicrobial finish.**

## 12 Year Warranty

5

*Our Manufacturer's Warranty is Unmatched* in the industry. We build our cabinets to last and we stand behind them.



## Superior Color Choices

10

*Over 20 Colors* to choose from. Cabinets have a heat sealed powder coat finish inside and out. View color options online at [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/)