

ACFF436L



81.5" x 54.25" x 29.75" (H x W x D)

39 cu.ft. general purpose healthcare freezer with French doors and stainless steel construction

Highlights:

- Made in North America in a TAA designated country
- Adjustable cantilevered shelving for flexible storage
- Audible/visual alarms for high/low temperature excursions

Product Features:

General purpose healthcare freezer made in North America	All-freezer designed to operate between -15 and -20.5°C for a variety of storage applications including laboratory specimens and freezer packs
TAA compliant	Ideal for GSA contracts and other applications requiring a unit to have been manufactured in a TAA designated country
Temperature Range: -15 to -20.5°C	Preset to -20°C as recommended for most items stored in general purpose applications
Frost-free	No-frost operation for minimum user maintenance
Fan-cooled interior	Powerful fan-forced cooling to ensure fast temperature recovery
High/low temperature alarm	Audible and visual alarm warns if the temperature goes out of range
Complete stainless steel construction	Interior and exterior are entirely constructed from stainless steel for lasting durability and easy cleanup
Self-closing doors	Doors automatically close if left ajar to help protect items in high traffic settings
Pocket handles	Handles recessed in the side of the doors for more efficiency in crowded facility settings

ACFF436L Specifications:

Overview	
Height of Cabinet	81.5" (207 cm)
Width	54.25" (138 cm)
Depth	29.75" (76 cm)
Depth with door at 90°	56.25" (143 cm)
Capacity	38.54 cu.ft. (1091 L)
Defrost Type	Frost-Free
Door	Stainless Steel
Cabinet	Stainless Steel
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Sanitation	UL-S
Amps	4.5
Voltage/Frequency	115 V AC/60 Hz
Weight	428.0 lbs. (194 kg)
Shipping Weight	450.0 lbs. (204 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101086101
Freezer	
Door Swing	French
Interior Height	55.0" (140 cm)
Interior Width	49.75" (126 cm)
Interior Depth	25.5" (65 cm)
Shelf Type	Wire
Shelf Qty	8
Adjustable Shelves	Yes
Thermostat Type	Digital
Interior Light	Yes
Fan Type	Interior
Refrigerant Type	R290a
Refrigerant Amount	4.2 oz.
High Side PSI	160.0
Low Side PSI	250.0
Wheels	6
Temperature Range	-20.5 to -15°C/-5 to +5°F

