

ACR1818



78" x 29.88" x 25.75" (H x W x D)

16.26 cu.ft. general purpose healthcare all-refrigerator with self-closing glass door, front lock, and digital controls

Highlights:

- Made in North America in a TAA designated country
- Reversible door swing to best accommodate your space
- Factory-installed lock included

Product Features:

General purpose healthcare refrigerator made in North America	All-refrigerator designed to operate between 2 and 8°C for a variety of storage applications including cold packs, laboratory specimens, and medications
Temperature Range: 2° and 8°C	Preset to 5°C as recommended for most items stored in general purpose applications
TAA compliant	Ideal for GSA contracts and other applications requiring a unit to have been manufactured in a TAA designated country
Glass door	Tempered glass door provides heat-safe view of interior contents
Self-closing door	Door automatically closes if left ajar to help protect items in high traffic settings
Factory-installed lock	Keyed lock for a secure interior (2 keys included)
Casters included	Pre-installed wheels allow you to easily move the unit as needed and lock into place with the front brakes
Pocket handle	Handle recessed in the side of the door for more efficiency in crowded facility settings
Fan-cooled interior	Powerful fan-forced cooling to ensure fast temperature recovery

ACR1818 Specifications:

Overview	
Height of Cabinet	78.0" (198 cm)
Width	29.88" (76 cm)
Depth	25.75" (65 cm)
Depth with door at 90°	55.25" (140 cm)
Capacity	16.26 cu.ft. (460 L)
Defrost Type	Automatic
Door	Glass
Cabinet	White
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Sanitation	UL-S
Amps	2.2
Voltage/Frequency	115 V AC/60 Hz
Weight	200.0 lbs. (91 kg)
Shipping Weight	250.0 lbs. (113 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101086088
Refrigerator Features	
Door Swing	RHD
Reversible	Yes
Shelf Type	Wire
Shelf Qty	3
Adjustable Shelves	Yes
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R290a
Refrigerant Amount	2.8oz.
High Side PSI	310.0
Low Side PSI	130.0
Wheel Qty	4
Interior Light	Yes
Temperature Range	2° to 8°C
Dimensions	
Interior Height	47.0" (119 cm)
Interior Width	25.25" (64 cm)
Interior Depth	21.88" (56 cm)

