

# S19LWH



20.13" x 17.75" x 17.63" (H x W x D)

Compact refrigerator-freezer with front-mounted lock for general purpose use; replaces S19L

### Highlights:

Designed for general purpose use in medical and scientific facilities

Factory-installed lock provides security you can depend on

Compact size is ideal for countertop use

### Product Features:

Factory installed lock	Keyed lock ensures a secure interior
Reversible door	User-reversible door swing allows for more flexibility in placing the unit
Manual defrost	Static cooling system
Door storage	Store large bottles right on the door for easy access
Adjustable thermostat	Dial thermostat conveniently located on the door to the freezer compartment for easy temperature management
Removable shelf	Wire shelf included for storage convenience
Compact dimensions	Sized for easy use on top of counters
100% CFC free	Environmentally friendly design without ozone-damaging chemicals
CARB Compliant	Eco-friendly hydrocarbon-based refrigerant for a greener footprint that meets current CARB regulations

## S19LWH Specifications:

Overview	
Height of Cabinet	19.75" (50 cm)
Height to Hinge Cap	20.13" (51 cm)
Width	17.75" (45 cm)
Depth	17.63" (45 cm)
Depth with door at 90°	33.88" (86 cm)
Capacity	1.7 cu.ft. (48 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	0.7
Voltage/Frequency	115 V AC/60 Hz
Weight	31.0 lbs. (14 kg)
Shipping Weight	56.0 lbs. (25 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101031910
Refrigerator-Freezer	
Door Swing	RHD
Reversible	Yes
Adjustable Shelves	Yes
Refrigerator - Interior Height	16.13" (41 cm)
Refrigerator - Interior Width	14.88" (38 cm)
Refrigerator - Interior Depth	14.25" (36 cm)
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	1
Refrigerator - Full Door Shelves	1
Refrigerator - Half Door Shelves	1
Freezer - Interior Height	4.38" (11 cm)
Freezer - Interior Width	8.75" (22 cm)
Thermostat Type	Dial
Refrigerant Type	R600a
Refrigerant Amount	0.42 oz.
High Side PSI	120.0

Low Side PSI	70.0
Compressor Step Height	7.38" (19 cm)
Compressor Step Width	14.88" (38 cm)
Compressor Step Depth	6.0" (15 cm)

