

FFAR23LGP



20.75" x 18.75" x 17.5" (H x W x D)

Compact 1.7 cu.ft. general purpose all-refrigerator in white with automatic defrost, front lock, and flat door liner

Highlights:

- Perfectly sized for countertop use in healthcare facilities
- Automatic defrost for easy user maintenance
- Factory-installed lock provides security you can depend on

Product Features:

Healthcare refrigerator	Designed for general purpose storage in institutional facilities
Countertop dimensions	Slim width and compact depth for easy placement on counters
Automatic defrost	Reduced maintenance with auto defrost system
Internal fan	Fan-forced cooling helps to keep the temperature more stabilized throughout the interior
Adjustable thermostat	Dial thermostat located inside unit for easy temperature management, with an operating range between 1 and 10°C
Front-mounted lock	Keyed lock for a secure interior (2 keys included)
Hospital grade cord with 'green dot' plug	Secure three-pronged cord for added unit safety in high traffic areas
LED lighting	Interior light for efficient illumination with an easy on/off rocker switch
Flat door liner	Allows the interior to store deeper items, while also providing improved temperature performance by eliminating hot spots commonly found in door shelves

FFAR23LGP Specifications:

Overview	
Height of Cabinet	20.25" (51 cm)
Height to Hinge Cap	20.75" (53 cm)
Width	18.75" (48 cm)
Depth	17.5" (44 cm)
Depth with door at 90°	35.0" (89 cm)
Capacity	1.7 cu.ft. (48 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.0
Voltage/Frequency	115 V AC/60 Hz
Weight	36.0 lbs. (16 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
Refrigerator Features	
Door Swing	RHD
Reversible	No
Shelf Type	Wire
Shelf Qty	1
Thermostat Type	Dial
Refrigerant Type	R600a
Refrigerant Amount	0.57oz.
High Side PSI	180.0
Low Side PSI	56.0
Level Legs Qty	2
Interior Light	Yes
Temperature Range	1 to 10°C
Dimensions	
Interior Height	16.25" (41 cm)
Interior Width	15.5" (39 cm)
Interior Depth	12.25" (31 cm)
Compressor Step Height	7.0" (18 cm)
Compressor Step Width	15.5" (39 cm)
Compressor Step Depth	6.5" (17 cm)

