

Product Name:	ABT-FFS-30
Description:	<p>General Purpose Flammable Storage Freezer 30 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Adjustable operating temperature range: -15°C – -25°C Four fixed shelves Keyed door lock Magnetic door gasket for positive seal Stainless steel door handle Leveling legs Probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation CFC-free refrigerant (R404a), foam insulation and packaging One year parts and labor warranty Exterior dimensions: 35 3/4"W x 35 1/2"D x 72 1/2"H Manual defrost 115V, 60Hz, 6 Amps, 1/3 HP UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 378 lbs.</p>