

# VT65ML



33.5" x 23.63" x 23.5" (H x W x D)

Freestanding medical all-freezer capable of -25° C operation with front-mounted lock

### Highlights:

-25°C operation for safe storage of temperature sensitive materials

Ideal 24" footprint for easy placement in institutional settings

Factory installed locks provide security you can count on

### Product Features:

-25°C operation	Low temperature operation is ideal for certain vaccines, pharmaceuticals, and other temperature sensitive materials
Factory installed lock	Keyed lock for a secure interior
Three removable storage drawers	Wire drawers with transparent plastic fronts allow for easy storage organization
Manual defrost	Reduced temperature fluctuations with an energy efficient manual defrost system
One-piece interior liner	Easy clean-up with a seamless liner that won't hold a mess in hidden crevices
Reversible door	User reversible door swing for added flexibility
Fully finished white cabinet	Allows the unit to be used freestanding
100% CFC Free	Environmentally friendly design without ozone-damaging chemicals

## VT65ML Specifications:

Overview	
Height of Cabinet	33.5" (85 cm)
Width	23.63" (60 cm)
Depth	23.5" (60 cm)
Capacity	3.5 cu.ft. (99 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Voltage/Frequency	115 V AC/60 Hz
Weight	105.0 lbs. (48 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101025100
Freezer	
Door Swing	RHD
Reversible	Yes
Interior Height	24.25" (62 cm)
Interior Width	17.75" (45 cm)
Interior Depth	16.5" (42 cm)
Thermostat Type	Dial
Refrigerant Type	R134a
Refrigerant Amount	2.76 oz.
Level Legs	2
Interior Drawers	Yes
Drawers Qty	3
Temperature Range	-12°C to -25°C
Sound Level	33 dB

