



Wood Style

MDS1830KC7-CM

M-SERIES TALL WOOD STYLE CART with Seven Drawers and a Key Lock

Why Wood Style Drawer Fronts?

Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.

Features:

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Sturdy, uni-body steel construction
- Vinyl coated aluminum drawer fronts
- Choice of lock options: www.harloff.com/locking-options/
- Full extension drawers slides for easy access
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Industry only fully curved protective bumper
- 30" of vertical drawer space with fully customizable drawer configurations
- Durable scratch resistant hard baked powder coat finish
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Pre-threaded mounting holes for accessories found at harloff.com/product-category/accessories/
- Latex free

Colors:

Shown with a Light Gray body, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

Wood Style drawer fronts are available in Cherry Mahogany or American Natural patterns, as shown to the right.



Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 42.625" (108.3 cm)
- Cart dims: 43.75" H x 24.75" W x 22" D (111.1 cm H x 62.9cm W x 55.9cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- Five 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada.
1-Year Limited Warranty on any bought-out parts.



Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of M-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff metal cart.



View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300