



# M-series

## MDSS3027B14 in Stainless Steel

### M-SERIES EMERGENCY CART Stainless Steel w/Four Drawers, Breakaway Lock

*Ideal for cleaning and sterilization, this cart includes the same great M-Series features but with a stainless steel body.\**

\* Not designed for use in a cart washer

#### Features:

- Sturdy, uni-body stainless steel exterior construction/finish with heat sealed powder coat steel drawer pans
- Choice of lock options
- Stainless steel breakaway lock bar
- Pack of 100 breakaway seals
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- 27" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters
- Latex free



Engineered & Manufactured in the U.S.



Shown in stainless steel. More than 20 additional colors are available in heat sealed powder coat steel.

#### Stainless Steel Medium Height, Standard Width, Four Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### Drawer Configuration:

- Two 3" drawers (7.6 cm), one 9" drawer (22.9 cm) and one 12" drawer (30.4 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

#### Assembled Cart Dimensions:

- 40.5" H x 36.75" W x 22" D (102.9 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 118 lbs (53.5 kg)
- Approximate shipping weight: 140 lbs (63.5 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at <https://www.harloff.com/acute-care-carts/emergency-crash-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300