



# M-series

## MDSS3027K16 in Stainless Steel

### M-SERIES TALL MEDICAL CART Stainless Steel w/Six Drawers and a Key Lock

*Ideal for cleaning and sterilization, this cart includes the same great M-Series features but with a stainless steel body.\**

\* Not designed for use in a cart washer

#### Features:

- Sturdy, uni-body stainless steel exterior construction/finish with heat sealed powder coat steel drawer pans
- Key lock with two keys
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- 27" of vertical drawer space
- Cart is delivered assembled
- Fully curved protective bumpers
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters

Shown in stainless steel. More than 20 additional colors are available in heat sealed powder coat steel.



#### Stainless Steel Tall, Standard Width, Six Drawer Medical Cart with Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



#### Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Three 3" drawers (7.6 cm) and three 6" drawers (15.2 cm)

**Engineered & Manufactured in the U.S.**

#### Assembled Cart Dimensions:

- 40.5" H x 36.75" W x 22" D (102.9 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 145 lbs (65.8 kg)
- Approximate shipping weight: 184 lbs (83.5 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at <https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300