

MPA1824E05+MD18-ANS



ALUMINUM

A-SERIES MINI WIDTH CART

w/Five Drawers and Electronic Pushbutton Lock



Features:

- Extremely light weight cart with critical grade aluminum construction, ideal for carts that are moved consistently throughout the day
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Antimicrobial plastic top, easy to remove/replace
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters
- Fully customizable drawer configurations at no extra charge
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at <https://www.harloff.com/locking-options/>
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 24" of vertical drawer space

Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.9 cm)
- 12.5" W x 17" D (31.8 cm x 43.2 cm)

Assembled Cart Dimensions:

- 37.25" H x 24.75" W x 22" D (94.6 cm H x 62.9 cm W x 55.9 cm D) without accessories attached
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm x 61 cm x 61 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Engineered & Manufactured in the U.S.



Color shown with a Light Gray body and Teal drawers. More than 20 additional colors are available.

Optional MD18-ANS Accessory Package (shown) Includes:

- 3-Tier Back Rail System (MD18-RBR3)
- (2) Aluminum Side Rails (SIDERAIL12)
- 5 Compartment Tilt Bin with Rail Clips (TILTBIN5RC)

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300