



MPA3030K16+MD30-AIRWAYPKG

ALUMINUM

A-SERIES TALL ANESTHESIA CART

6 Drawer Key Lock, MD30-AIRWAYPKG Anesthesia Pkg

MD30-AIRWAYPKG Accessories Include:



Oxygen tank holder (O2HLDR)



(2) Utility hooks (UHOOKS2)



(1) 16" Side Rail (SIDERAIL16)



Adjustable Overhead Shelf (MD30-ADJSHLF2)



5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)



4-Scope Capacity Side Scope Cabinet (680350)

Drawer Configuration:

- Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.4 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

Assembled Cart Dimensions:

- 66.75" H x 41.25" W x 22" D (169.5 cm H x 104.78 cm W x 55.9 cm D)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Ergonomic side push handles, mounted on the left or right
- Antimicrobial plastic top, removable and easily replaceable
- Fully customizable drawer configurations at no extra charge
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 30" of vertical drawer space

Color shown is Navy. More than 20 additional colors are available.



View similar products online at <https://www.harloff.com/acute-care-carts/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com



Engineered & Manufactured in the U.S.