

Stainless Steel Cast Cart/ Model: 6020E

Eight Drawer with Basic Electronic Push Button Lock, Standard Package



Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

Construction:

- Stainless steel cabinet and drawers

Features:

- Stainless steel bucket holder with vinyl mat
- 24" stainless steel ruler etched on top of cart
- Aluminum towel bar
- Internal drawer dividers
- Protective corner bumpers
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Latex free
- Designed in conjunction with orthopedic technicians to maximize functionality

Other Cast Cart options:

- All cast carts are also available with three different lock types: breakaway, key-locking, and basic electronic pushbutton locks. Call Harloff for details.

Drawer Configuration:

- (1) - 3-inch with black vinyl mat (7.6 cm)
- (5) - 3-inch with 4 adj. dividers per drawer (7.6 cm)
- (2) - 6-inch with 4 adj. dividers per drawer (22.9 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Overall Cart Dimensions:

- 43" H x 38" W x 21" D (109.2 cm H x 96.5 cm W x 53.3 cm D)
- Writing surface: 42" H x (106.7 cm)
- Cart net weight: 152 lbs (69.1 kg)
- Approximate shipping weight: 240 lbs (109.1 kg)
- Approximate shipping dims: 57" H x 47" W x 27" D (144.8 cm H x 119.4 cm W x 68.6 cm D)