



R.

M-SERIES PPE & ISOLATION CART

Four Drawers with a Standard Key Lock

Features:

- Includes a sani-wipe holder accessory (SANIHLDR) and a 3-box glove holder accessory (GLOVEBOX3DM)
- · Space-saver design fits into smaller areas allowing supply storage in areas with space limitations
- · Sturdy, uni-body steel construction
- · Key lock with two keys
- · Durable scratch resistant heat sealed powder coat finish
- · Full extension drawers slides for easy access
- · Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 30" of vertical drawer space
- · Cart is delivered assembled
- · Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free

Warranty:

· 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.

Color shown in Yellow. More than 20 additional colors are available.

Gloves and sani-wipes must be provided by your facility.

Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 39.625" (100.6 cm)
- Cart dims: 40.75" H x 18" W x 18" D (103.5 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

One 3" drawer (7.6 cm) and three 9" drawers (22.8 cm)

Designed

and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



View similar products online at https://www.harloff.com/product-category/m-series-medical-carts/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300