



# M-series

## MDS1824K05+MD18-PHB

### M-SERIES SHORT PHLEBOTOMY CART Mini Five Drawer Cart w/MD18-PHB Accessory Package

#### Short Height, Mini Width, Five Drawer Phlebotomy Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Phlebotomy carts keep supplies nearby, are mobile, provide a dedicated work surface for procedures allow quick disposal of used needles. M-Series carts offer you choices.

#### MD18-PHB Accessories Include:



6-Compartment Tilt Bin w/Rail Clips (TILTBIN6RC)



Glove Box Holder Single Direct Mount (GLOVEBOX1DM)



Aluminum Waste Container with Lid (AL680436)



Push Handle (MINI-PBHANDKIT)



3-Tier Back Rail (MD18-RBR3)



Engineered & Manufactured in the U.S.

Color shown in Hammertone Blue. More than 20 additional colors are available.



- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations within 24" vertical space
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters



View similar products online at

<https://www.harloff.com/product-category/phlebotomy-carts-iv-start-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300  
medicalsales@harloff.com • www.harloff.com

#### Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Cart dims: 36.38" H x 23.13" W x 22" D (92.4 cm H x 58.7 cm W x 55.9cm D)
- Cart net weight: 86 lbs (39 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm H x 61 cm W x 61 cm D)

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) & one 9" drawer (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.