



M-series

MDS1830K06+MD18-PHB

M-SERIES TALL PHLEBOTOMY CART Mini Six Drawer Cart w/MD18-PHB Accessory Package

Tall Height, Mini Width, Six Drawer Phlebotomy Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Phlebotomy carts keep your supplies nearby, are mobile, provide a dedicated work surface for procedures and allow quick disposal of used needles. M-Series carts offer you choices.

MD18-PHB Accessories Include:



6-Compartment Tilt Bin w/Rail Clips (TILTBIN6RC)



Glove Box Holder Single Direct Mount (GLOVEBOX1DM)



Aluminum Waste Container with Lid (AL680436)

Engineered & Manufactured in the U.S.



Push Handle (MINI-PBHANDKIT)



3-Tier Back Rail (MD18-RBR3)

Color shown with a Cream body and White drawers. More than 20 additional colors are available.



- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawer slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Cart is delivered assembled
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters
- Fully customizable drawer configs with 30" of vertical drawer space



View similar products online at

<https://www.harloff.com/product-category/phlebotomy-carts-iv-start-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Overall Cart Dimensions:

- Work Surface Height: 41.5" (105.4 cm)
- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Cart dims: 42.625" H x 23.875" W x 22" D (108.3 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.