



M-series

MDS1830L06+MD18-PHB

M-SERIES TALL PHLEBOTOMY CART Mini Six Drawer Cart w/MD18-PHB Accessory Package

Tall Height, Mini Width, Six Drawer Phlebotomy Cart with a Latch Lock

Designed and built to meet the needs of today's healthcare market. Phlebotomy carts keep your supplies nearby, are mobile, provide a dedicated work surface for procedures and allow quick disposal of used needles. M-Series carts offer you choices.

MD18-PHB Accessories Include:



6-Compartment Tilt Bin w/Rail Clips (TILTBIN6RC)



Glove Box Holder Single Direct Mount (GLOVEBOX1DM)



Aluminum Waste Container with Lid (AL680436)

Engineered & Manufactured in the U.S.



Push Handle (MINI-PBHANDKIT)



3-Tier Back Rail (MD18-RBR3)

Color shown in Navy. More than 20 additional colors are available.



- Sturdy, uni-body steel construction
- Turn latch lock to secure drawers while moving cart
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawer slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Cart is delivered assembled
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters
- Fully customizable drawer configs with 30" of vertical drawer space



View similar products online at

<https://www.harloff.com/product-category/phlebotomy-carts-iv-start-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Overall Cart Dimensions:

- Work Surface Height: 41.5" (105.4 cm)
- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Cart dims: 42.625" H x 23.875" W x 22" D (108.3 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.