



# MDS2430KC4+MD24-DRW18-EHDR

*Equipment Holder with Door*

## M-SERIES MEDIUM WIDTH CART

2 Drawer, 18" Equipment Holder w/Door & Key Lock

### M-Series Features:

- Sturdy, uni-body steel construction
- Choice of lock options:  
[www.harloff.com/locking-options/](http://www.harloff.com/locking-options/)
- Up to 30" of fully customizable vertical space is available for both drawer and equipment holder configurations
- 17" deep closed door equipment holding area can replace standard drawer areas in three possible height options
- 2" pass-through hole in the rear for cords and cables
- Full extension drawer slides for easy access
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Antimicrobial plastic top, removable
- Durable scratch resistant heat sealed powder coat finish
- Body has pre-threaded accessory mounting holes
- Industry best 5" quiet rolling Tente® casters  
- 1 steering caster and 2 braking casters
- Latex free



**MD24-DRW18-EHDR**  
**Equipment Holder**

**Color shown is Teal.**  
**More than 20 additional colors are available.**



**Engineered & Manufactured in the U.S.**

**M-series**

#### Overall Cart Dimensions:

- Internal Drawer Dimensions: 18.5" W x 17" D (47 x 43.2 cm)
- Work Surface Height: 42.625" (108.3 cm)
- Cart dims: 43.75" H x 30" W x 22" D (111.1 cm H x 76.2cm W x 55.9cm D)
- Cart net weight: 119 pounds (54 kg)
- Approx shipping weight: 155 pounds (70.3 kg)
- Approx shipping dims: 56" H x 33" W x 26" D (142.2 cm x 99 cm x 66 cm)

#### Drawer Configuration:

- Two 6" drawers (15.2 cm) and one MD24-DRW18-EHDR doored equipment holder with an 17.4" H x 18.63" W x 17" D (44.1 cm H x 47.3cm W x 43.2cm D) opening

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Designed and built to meet the needs of today's healthcare market. Replace 9", 12 or 18" of vertical drawer space at the top or bottom of your M-Series cart with a doored equipment holder, perfect for storing electronics or large medical equipment. M-Series carts offer you choices.

View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)**  
**[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**

**HARLOFF**  
INNOVATING STORAGE FOR CARE