



MDS3024SC

BABY SCALE CART

with One Drawer Plus an Open Holder with Shelf

The cart top is designed to store a baby scale, plus one drawer and an adjustable shelf for extra equipment and supply storage

Supplies not included

Features:

- Sturdy, uni-body steel construction
- Ideal for holding a baby scale on the cart top and other large equipment in the interior
- Includes one 3" drawer plus 21" of open vertical shelf space
- Includes adjustable shelf with seven possible height locations
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawer slides for easy access
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Body has pre-threaded accessory mounting holes
- Ergonomic side push handles can be mounted on either side
- Cart is delivered assembled
- Latex free

Color shown in Beige. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

Short Height, Standard Width Baby Scale Cart with One Drawer and an Adjustable Shelf

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls and expanded lock options make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

Drawer Configuration:

- One 3" drawer (15.2 cm) plus an open equipment holder with an 22" H x 25.75" W x 17" D (51.8 cm H x 62.6cm W x 43.2cm D) opening with adjustable shelf
- 24.5" W x 17" D (62.2 x 43.2 cm)

Assembled Cart Dimensions:

- 35.89" H x 36" W x 22" D (91.1 cm H x 91.4 cm W x 55.9 cm D)
- Cart weight: 71 lbs (32.2 kg)
- Approximate shipping weight: 110 lbs (49.9 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts.



View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300