



# M-series MDS3030B09PED

## M-SERIES TALL PEDIATRIC CART with Nine Drawers & Breakaway Drawer Locks

### Features:

- Sturdy, uni-body steel construction
- Tall cabinet with separate breakaway locks on each drawer
- Cart is designed to be compatible with popular pediatric color coding systems
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Ergonomic side push handles can be mounted on either side
- Comfort drawer pulls
- Industry best 5" quiet rolling Tente® casters - 1 steering caster & 2 braking casters
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Painted steel break-away lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Latex free
- Cart is delivered assembled



Engineered &  
Manufactured  
in the U.S.

**Color shown is Beige with color coded drawers. More than 20 additional body colors are available.**



### Tall, Standard Width, Nine Drawer Pediatric Cart

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### Drawer Configuration:

- Eight 3" (7.6 cm) drawers and one 6" (15.2 cm) drawer
- 24.5" W x 17" D (62.2 x 43.2 cm)

#### Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (1111 cm H x 93.3 cm W x 55.9 cm D)
- Writing surface: 42.625" (108.2 cm)
- Cart weight: 144 lbs (65.3 kg)
- Approximate shipping weight: 185 lbs (83.9 kg)
- Approximate shipping dims: 56" H x 39" W x 26" D (142.2 cm x 99 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at

<https://www.harloff.com/product-category/pediatric-resuscitation-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300