



M-series MDS3030B09PED1

M-SERIES TALL PEDIATRIC CART with Nine Drawers and a Breakaway Lock

Designed for treating pediatric emergencies by using a multi-colored drawer labeling system comparable to the Broselow color coding method

Features:

- Sturdy, uni-body steel construction
- Heat sealed powder coat steel full length breakaway lock bar
- Pack of 100 breakaway seals
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Ergonomic side push handles can be mounted on either side
- Comfort drawer pulls
- Industry best 5" quiet rolling Tente® casters - 1 steering caster & 2 braking casters
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 30" of vertical drawer space
- Latex free
- Cart is delivered assembled

Color shown is Light Gray with color coded drawers. More than 20 additional body colors are available.



Tall, Standard Width, Nine Drawer Pediatric Cart

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

Drawer Configuration:

- Eight 3" (7.6 cm) drawers and one 6" (15.2 cm) drawer
- 24.5" W x 17" D (62.2 x 43.2 cm)

Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (1111 cm H x 93.3 cm W x 55.9 cm D)
- Writing surface: 42.625" (108.2 cm)
- Cart weight: 144 lbs (65.3 kg)
- Approximate shipping weight: 194 lbs (88 kg)
- Approximate shipping dims: 56" H x 39" W x 26" D (142.2 cm x 99 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

**Engineered &
Manufactured
in the U.S.**

View similar products online at

<https://www.harloff.com/product-category/pediatric-resuscitation-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**

HARLOFF
INNOVATING STORAGE FOR CARE