



M-series

MDSS2424B13 in Stainless Steel

M-SERIES SHORT EMERGENCY CART with Three Drawers and a Breakaway Lock

*Ideal for cleaning and sterilization, this cart includes the same great M-Series features but with a stainless steel body.**

* Not designed for use in a cart washer

Features:

- Sturdy, uni-body stainless steel exterior construction/finish with heat sealed powder coat steel drawer pans
- Stainless steel breakaway lock bar (other locking options are available)
- Pack of 100 breakaway seals
- 24" of vertical drawer space (3", 9" and 12" high drawers)
- Full extension drawer slides for easy access
- Body has pre-threaded accessory mounting holes
- Ergonomic side push handles can be mounted on either side
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - one steering caster and two braking casters
- Latex free

Shown in stainless steel. More than 20 additional colors are available in heat sealed powder coat steel.



Short Height, Medium Width, Three Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered & Manufactured in the U.S.

Drawer Configuration:

- One 3" drawers (7.6 cm), one 9" drawer (22.8 cm) and one 12" drawer (30.4 cm)
- 18.5" W x 17" D (47 cm x 43.2 cm)

Assembled Cart Dimensions:

- 37.25" H x 30" W x 22" D (94.6 cm H x 76.2 cm W x 55.9 cm D)
- Cart net weight: 111 pounds (50.3 kg)
- Approx shipping weight: 150 pounds (68 kg)
- Approx shipping dims: 50" H x 33" W x 26" D (127 cm x 99 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300