



Key Locking

DSC24SK

SHORT DILATOR STORAGE CART for Savary-Galliard® Dilators up to 28" in length

Current guidelines require HEPA filtration and require that endoscopes be hung vertically due to the need for inner channels to drain all excess fluid after cleaning. Multiple facilities have requested an esophageal dilator storage cabinet to allow them to store their Savary-Gilliard® Dilators vertically because they also have an internal channel.

With this in mind, Harloff has introduced the new Dilator Storage Cart. Vertically hang up to 24 dilators for gravity drying and Savary dilator storage. Our standard esophageal dilator storage cabinet can hold dilators up to 28" in length.



Engineered & Manufactured in the U.S.

Color shown is White with a Teal Door.
More than 20 additional colors are available.



Features:

- Each Savary dilator storage cart can hold up to 24 individual dilators, varying in size from 5mm (15 French) to 20mm (60 French) and up to 28" in length
- Includes 4 Savary pull out slides that each hold six dilators for easy access to similar sized dilators at one time
- Choice of lock options: latch lock, standard key lock (shown) or standard electronic lock
- Sturdy, fully welded steel body
- Industry best 5" quiet rolling Tente® casters - 1 steering caster & 2 braking casters
- Industry only fully curved protective bumper
- Durable scratch resistant heat sealed powder coat finish

Assembled Cart Dimensions:

- 40.38" H x 23.89" W x 22" D (102.6 cm H x 60.7 cm W x 55.9 cm D)
- Cart weight: 85 lbs (38.6 kg)
- Approximate shipping weight: 110 lbs (49.9 kg)
- Approximate Shipping Dims: 50" H x 29" W x 24" D (127 cm x 73.7 cm x 61 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts



HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/scope-drying-and-storage-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300