



Key Locking

DSC24TK

TALL DILATOR STORAGE CART for Savary-Gilliard® Dilators up to 39" in length

Another new medical innovation from the design and engineering team at Harloff!

Achieve the same storage capabilities of an endoscope cabinet with your Savary-Gilliard® Dilators. This new vertical dilator storage cabinet can accommodate up to 24 dilators for storage and gravity drying. Our tall Savary dilator storage cart can hang dilators up to 39" in length, and comes with your choice of four lock options.



Engineered & Manufactured in the U.S.



Color shown is Light Gray with a Navy Door. More than 20 additional colors are available.

Features:

- Each Savary dilator storage cabinet can hold up to 24 individual dilators, varying in size from 5mm (15 French) to 20mm (60 French) and up to 39" in length
- Includes 4 Savary pull out slides that each hold six dilators for easy access to similar sized dilators at one time
- Choice of lock options: latch lock, standard key lock (shown) or standard electronic lock
- Sturdy, fully welded steel body
- Industry best 5" quiet rolling Tente® casters - 1 steering caster & 2 braking casters
- Industry only fully curved protective bumper
- Durable scratch resistant heat sealed powder coat finish



Assembled Cart Dimensions:

- 51.89" H x 23.89" W x 22" D (131.8 cm H x 60.7 cm W x 55.9 cm D)
- Cart weight: 97 lbs (44 kg)
- Approximate shipping weight: 122 lbs (55.3 kg)
- Approximate Shipping Dims: 50" H x 29" W x 24" D (127 cm x 73.7 cm x 61 cm)

Warranty:

12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts



View similar products online at <https://www.harloff.com/scope-drying-and-storage-cabinets/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300