

# ARFG18PVDL2B



79.63" x 27.63" x 26.25" (H x W x D)

PHARMA-VAC Performance Series combination 2-door/2-compressor refrigerator-freezer with factory-installed digital data loggers for pharmacy, medical, and life science applications

## Highlights:

Purpose-built design to store medication, pharmaceuticals, and other sensitive materials

Independent dual compressor refrigerator and freezer storage in one slim-fitting footprint to maximize storage capacity without compromising performance

Intelligent microprocessor digital temperature controller for optimum temperature management

## Product Features:

Performance Series Refrigerator-Freezer Combination	Purpose-built all-refrigerator and all-freezer for vaccine storage in medical, pharmaceutical, and life science applications
Purpose-built design	Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC/VFA guidelines
+/-1°C Variation	The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period
Digital data loggers (DDL) included	Accucold DL2B data logger is NIST calibrated in our ISO/IEC 17025:2017 laboratory, with user-set logging intervals, a large LCD readout of the current and high/low temperature, and virtually unlimited storage through internal memory and a convenient USB drive; for additional information on this logger's features, view the DL2BKit on our website
Refrigerator window	Triple paned glass on the upper refrigerator section provides a clear view of interior contents
Space-saving design	Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint
Antimicrobial silver ion-coated handles	Powder-coated handles with naturally occurring silver ions to help reduce the spread of germs
Adjustable temperature	Intuitive controls let you adjust the interior temperature between +2 to +8°C (refrigerator) and -30 to -20°C (freezer)



## ARFG18PVDL2B Specifications:

Overview	
Height of Cabinet	79.63" (202 cm)
Width	27.63" (70 cm)
Width with Door Open	28.13" (71 cm)
Depth	26.25" (67 cm)
Depth with Handle	27.25" (69 cm)
Depth with door at 90°	55.13" (140 cm)
Capacity	15.54 cu.ft. (440 L)
Defrost Type	Cycle
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	2.6
Voltage/Frequency	115 V AC/60 Hz
Weight	275.0 lbs. (125 kg)
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
UPC	761101159744
Refrigerator-Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	40.0" (102 cm)
Refrigerator - Interior Width	23.75" (60 cm)
Refrigerator - Interior Depth	17.75" (45 cm)
Refrigerator - Capacity	11.04 cu.ft
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	4
Refrigerator - Temperature Range	2 to 8°C
Refrigerator - Refrigerant Type	R600a
Refrigerator - Refrigerant Amount	1.44 oz.
Freezer - Interior Height	15.75" (40 cm)
Freezer - Interior Width	22.0" (56 cm)

Freezer - Interior Depth	19.5" (50 cm)
Freezer - Capacity	4.5 cu.ft
Freezer - Shelf Type	Wire
Freezer - Shelf Qty	2
Thermostat Type	Digital
Fan Type	Interior
Wheels	4
Freezer - Temperature Range	-30 to -10°C
Freezer - Refrigerant Type	R290a
Freezer - Refrigerant Amount	1.27 oz.

